

# CuroCell Cirrus<sup>®</sup> 2.0



Pain therapy pressure care system

# CuroCell Cirrus® 2.0

CuroCell Cirrus® 2.0 is used as an aid in the prevention and treatment of pressure ulcers up to category 4 (0-180 kg) and category 2 (181-230 kg)<sup>(1)</sup>, and for pain therapy. It is a replacement system ideal for all healthcare environments.

CuroCell Cirrus® 2.0 combines air with a soft foam surface and provides high user comfort and pressure redistribution<sup>(2)</sup>. The unique wave like alternating process is particularly suitable for pain therapy. The added heel function offers an even better protection by relieving pressure on the heels.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene<sup>(3)</sup>, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(4)</sup>, which reduce the risk of skin maceration.

## Facts

Type: Alternating pressure mattress

Height: 21 cm

Function: Prevention and treatment up to and including category 4<sup>(1)</sup>

Recommended user weight: 0-180 kg (cat. 4) & 180-230 kg (cat. 2)

Fire: EN 597-1, -2

Power supply: 230 V/50 Hz

Noise level of pump: 27 dBA

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



## Features and Characteristics

- Alternating pressure on soft foam surface
- Digital and easy to use pump
- Easy handling and cleaning
- Heel slope, for pressure relief of heels
- Hygienic cover with welded seams
- Rapid deflation of the mattress
- Silent and energy-saving pump
- Soft, elastic, vapour permeable<sup>(4)</sup> hygienic cover
- Sound and light alarm for low pressure
- Sound and light alarm for power failure
- Static function with automatic return
- Transport function

The silent and easy to use CuroCell Cirrus® pump with four different settings; alternating, static, alternating + seat inflation and static + seat inflation make this mattress very functional for all patient groups.

In addition the dual cover system reduces the shearing effect when adjusting the bed.



## Sizes

80/85/90/105/120x200 cm and 80/85/90x210 cm.

If you need any other sizes please contact your distributor.

## Optional

- Durable bag
- Low friction cover

## References

(1) European Pressure Ulcer Advisory Panel, Quick reference guide 2009.

(2) Examination of the shear force and the pressure relief capability, Berlincert 2009.

(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(4) SS-EN ISO 15496.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell Cirrus are trademarks of Care of Sweden AB.

© Care of Sweden AB, 2014. All rights reserved. Last update: 26/08/14.



Scan QR code to watch presentation video.



Scan QR code to watch instructional video.



### Mattress selector

- Guidelines when choosing a product  
See how to choose the right mattress at:  
[www.careofsweden.com/mattress-selector](http://www.careofsweden.com/mattress-selector)

 **care of sweden**  
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden  
+46 (0)771 106 600 | [info@careofsweden.se](mailto:info@careofsweden.se) | [www.careofsweden.com](http://www.careofsweden.com)