

CuroCell UNO®



A new generation of fully automated
alternating pressure mattress replacement systems

CuroCell UNO®

The CuroCell UNO® is used as an aid in the prevention and treatment of pressure ulcers up to and including category 4⁽¹⁾. It stands for innovation and modern design and together with quality engineering it lays the foundation for our new generation of fully automated premium mattress systems.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽³⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which reduce the risk of skin maceration.

Facts

Type: Alternating pressure mattress replacement system
Height: 20 cm
Function: Prevention and treatment up to and including category 4⁽¹⁾
Recommended user weight: 0-250 kg
Input voltage: 100-240 V / 50-60 Hz
Sound level, pump, max.: 20 dBA⁽²⁾
Warranty: 2 years
Cleaning of cover: Wipe with cleaning agent and/or disinfectants
Machine wash max 95 °C, tumble drying



Features Pump

- Fully automated pump that adapts individually
- Large colour touch display with user-friendly interface
- 3 modes: Alternating, Pulsating and Static
- Pack&Go™ function
- Silent running and vibration free
- Battery backup (optional)
- Soft friction hangers with extra soft tip
- Notification function for loss of air and or electricity supply

Features Mattress

- Cell on Cell construction with safety mattress
- Heel Function
- Static head cells
- CPR-function for rapid deflation
- Smart cable management
- PVC free
- Unique design for easy handling and usage
- Chlorine resistant (10 %) multi stretch cover with welded seams
- Well protected zipper

It is equipped with a large touch screen panel with a user-friendly interface. The fully automated pump offers automatic adjustment, providing optimal pressure relief and comfort for each individual user while aiding simple use.

Sizes

80/85/90/105/120 x 200/210/220 x 20 cm
If you need any other sizes please contact your distributor

Optional

- Durable bag
- Sheet holder

References

- (1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.
(2): "EN ISO 11201 Acoustics - Noise emitted by machinery and equipment - Determination of emission sound pressure levels", SP 2015.
(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles - Bacteria penetration - Wet".
(4) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell UNO are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2016. All rights reserved. Last update: 2016-01-19.



Scan QR code to watch presentation video.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 **care of sweden**
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com