

CuroCell S.A.M.® PRO

Extraordinary comfort
and durability combined



CuroCell S.A.M.[®] PRO

The CuroCell S.A.M.[®] PRO overlay mattress can be used as an effective aid in the prevention and in the treatment of pressure ulcers up to and including category 3⁽¹⁾. It can also be used as an aid for pain therapy. By using the optional inner cover comfort can be increased while further reducing pressure.

This static air mattress offers high level of comfort as well as effective pressure⁽²⁾ redistribution. With the vastly improved higher maximum weight capacity of 250 kg, this overlay will meet the ever increasing challenges of heavier patients.

The CuroCell S.A.M.[®] PRO features a new cost effective pump that includes both CPR and Transport functions. The pump automatically determines when to change pressure in contrast to running constantly. This makes the pump silent running, vibration free, power efficient and long lasting.


Also available is the CuroCell S.A.M.[®] PRO CF16, which is a mattress replacement system equipped with an integrated foam core, which is placed directly onto the bed frame.



Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽³⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which reduce the risk of skin maceration.

Facts

Type: Static air mattress
Height: 10 cm
Function: Prevention and treatment up to and including category 3⁽¹⁾
Recommended user weight: 0-250 kg
Warranty: 2 years
Cleaning of cover: Wipe with cleaning agent and/or disinfectants.
Machine wash max 95 °C, tumble drying.



Features & Characteristics

- CPR function – rapid deflation
- Easy to handle and maintain
- Elastic fastening band
- Silent running and durable pump
- Energy efficient
- Static function and Max firm function
- Comfort settings on dial, 0-250 kg
- Soft, elastic, vapour permeable⁽⁴⁾ hygienic cover with welded seams
- Unique labelling
- Very good cleaning properties

Sizes

80/85/90/105/120x200/210 cm.
If you need any other sizes please contact your distributor.

Optional

- Durable bag for transportation and storage
- Extra soft inner cover
- Transport function

References

- (1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.
(2) Examination of the shear force and the pressure relief capability, Berlincert 2009.
(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".
(4) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. are trademarks of Care of Sweden AB. © Care of Sweden AB, 2017. All rights reserved. Last update: 2017-04-18.



Scan QR code to watch presentation video.



Mattress selector

- Guidelines when choosing a product

See how to choose the right mattress at:

www.careofsweden.com/mattress-selector

 **care of sweden**
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com